

OMRON

Digital Body Weight Scale

HN-283 MODEL



Measures body weight, in 100g increment. Even a slight weight change can be detected.

Easy operation, just tap to switch on

Sleek and trendy design

Silver coated advanced tempered glass

4 sensors for greater precision

Large LCD panel

Specifications

Model	Digital Body Weight Scale HN-283	
Display	Body weight: 5-150kg with an increment of 100g	
Weight accuracy	5 kg - 60 kg: $\pm 400g$ / 60.1 - 100kg: $\pm 600g$ / 100.1kg - 150kg: $\pm 1\%$	
Power	Two 3V lithium batteries	
Battery life	Approximately 1 year (3 times measurements per day)	
Weight	Approximately 1.4kg (excluding batteries)	
External dimensions	W x H x D (mm)	300(W) x 320(D) x 40(H)
Accessories	2 test batteries, User Guide	

NOTE: These specifications are subject to change without notice. Printed 2008

Distributed by:

A Good Sense of Health

OMRON HEALTHCARE SINGAPORE PTE LTD
83 Clemenceau Avenue, #15-05, UE Square, Singapore 239920
Website: www.omron-healthcare.com.sg